



# 2024

*Product Catalogue*

# Company Profile

## About Us

Central Medical (Pty) Ltd was established in 1996 as a local manufacturer of disposable non-woven products, arguably one of the first to introduce automated manufacturing processes for medical supplies into South Africa.

The Company has since diversified its operations to also perform direct importation of additional hygienic disposable products. Over the years the Company has developed a number of internal brands that has accumulated significant reputation for being reliable, consistent and affordable. These brands include LYSA, Centromed, Hygienic & Protective (H&P) and Examtex.

## What we can help you with

In addition to providing goods currently listed in our product catalogue, we are able to provide you with medical supplies in your specified branding as well as any medical supplies that would further suite your needs.

## What we stand for

### *Comprehensive understanding of our offerings*


Central Medical has manufactured latex gloves, non-woven apparels and bandages locally. In addition, we have long term relationships with our overseas suppliers since the company's inception. The experience held by our team has allowed us to ensure that the products manufactured would meet manufacturing and regulatory standards.

### *Consistent offerings*

We uphold the same manufacturing standards which we implement locally to our overseas long-term suppliers in order to ensure that the quality of our products remain consistent.

### *Affordable and Competitive*

As a manufacturer and a direct importer, Central Medical avoids the additional costs charged by agents. This cost benefit is reflected in our prices, meaning that it is passed onto our partners.

  
CEO



# Contact

*Central Medical Branches*

# Branches



## Durban – Head Office

28 Duiker Road, Canelands, 4339

[sales@centromed.co.za](mailto:sales@centromed.co.za)

T: +27 32 533 9724

C: +27 82 800 6666



## Johannesburg – Branch

88 Cranberry Crescent, Unit 6 Acacia Park

[sales-jhb@centromed.co.za](mailto:sales-jhb@centromed.co.za)

C: +27 73 493 2625

C: +27 64 065 3470



A photograph of industrial machinery, likely a paper mill, featuring a series of large rollers and mechanical components. The image is partially obscured by a blue overlay and a white box.

2

# Products

*2024 Offerings*



# Needles & Syringes

*Sterile, Single-use only*

# Hypodermic Needle

*Sterilised, Single-use*



## Product Features

- Thin-walled needles
- Made of austenite stainless steel
- Smooth surface with light lubricative coating
- Minimal pain upon puncture
- Transparent Luer-lock plastic hub
- Color-coded, medical grade plastic hubs for easy identification
- E.O. Sterilised
- Non-pyrogenic
- Box of 100 pieces

| Code No.  | Gauge x Length | Ø x Length (mm) | Colour Coding | Units /Carton |
|-----------|----------------|-----------------|---------------|---------------|
| NED1HPD01 | 15G x 1,5"     | 1,8 x 40        | Blue-grey     | 100           |
| NED1HPD02 | 18G x 1,5"     | 1,2 x 40        | Pink          |               |
| NED1HPD03 | 19G x 1,5"     | 1,1 x 40        | Cream         |               |
| NED1HPD04 | 20G x 1,5"     | 0,9 x 40        | Yellow        |               |
| NED1HPD05 | 21G x 1,5"     | 0,8 x 40        | Deep-green    |               |
| NED1HPD06 | 22G x 1,25"    | 0,7 x 32        | Black         |               |
| NED1HPD07 | 23G x 1"       | 0,6 x 25        | Deep-blue     |               |
| NED1HPD09 | 25G x 0,625"   | 0,5 x 16        | Orange        |               |
| NED1HPD10 | 26G x 0,625"   | 0,45 x 16       | Brown         |               |
| NED1HPD11 | 29G x 0,33"    | 0,33 x 8        | Red           |               |



# Needles

*Sterilised, Single-use*





## Product Features - Blood Lancets

- Twist-off protective cap
- Pressure activated
- Gamma Sterilised
- Single-use only
- 100 pieces per pack

| Code No.  | Gauge       | Colour Coding  | Units /Carton |
|-----------|-------------|--|---------------|
| LAN1ECO04 | 21G * 1,8mm | Green   | 20            |
| LAN1ECO02 | 25G * 1,6mm | Yellow  |               |
| LAN1ECO03 | 28G * 1,4mm | Orange  |               |

## Product Features – ReadyLance Blood Lancets

- Twist-off protective cap
- Pressure activated
- Gamma Sterilised
- Single-use only
- 100 pieces per pack

| Code No.  | Gauge       | Colour Coding  | Units /Carton |
|-----------|-------------|--|---------------|
| LAN1STD01 | 25G * 1,6mm | Yellow  | 20            |
| LAN1STD02 | 28G * 1,4mm | Orange  |               |







# Needles

*Sterilised, Single-use*



## Product Features – TB Syringes

- Thin-walled long-point needle
- Highly transparent medical polypropylene (PP) barrel
- E.O. Sterilised
- Non-pyrogenic
- Manufactured to ISO 8537 requirements
- 100 pieces per pack

| Code No.  | Capacity | Gauge      | Packs /Carton |
|-----------|----------|------------|---------------|
| SYR1TBA01 | 1 ml     | 25G x 5/8" | 20            |





# Insulin syringes

*Sterilised, Single-use*



## Product Features

- Double capped to ensure sterility with ultra-thin needles
- No dead-space within syringe cavity to reduce insulin wastage
- Ultra thin, sharp and lubricated needle.
- E.O. sterilised
- Non-pyrogenic
- U-100/U-50
- Manufactured to ISO 8537 requirements
- 100 pieces per pack

| Code No.  | Capacity | Gauge    | Packs /Carton |
|-----------|----------|----------|---------------|
| SYR1INS01 | 0,3 ml   | 29G x ¼" | 20            |
| SYR1INS02 | 0,5 ml   | 29G x ¼" |               |
| SYR1INS03 | 1,0 ml   | 29G x ¼" |               |
| SYR1INS04 | 0,5 ml   | 27G x ½" |               |
| SYR1INS05 | 1,0 ml   | 27G x ½" |               |
| SYR1INS06 | 1,0 ml   | 29G x ½" |               |
| SYR1INS07 | 1,0 ml   | 26G x ½" |               |








# Syringes

*Sterilised, Single-use, Luer-Lock*



## Product Features

- Transparent medical polypropylene (PP) barrel.
- Luer-lock tip
- Automated assembly process to prevent personnel contamination.
- Synthetic rubber stopper to reduce allergic reactions.
- Medical grade silicone lubrication to improve sealing ability and glide.
- E.O. Sterilised
- Non-pyrogenic

| Code No.  | Capacity   | Concentr. /Eccentr. | Colour Coding |   | Units /Box |
|-----------|------------|---------------------|---------------|---|------------|
| SYR1HPL01 | 2 / 3 ml   | Concentric          | Blue          |  | 100        |
| SYR1HPL02 | 5 ml       | Concentric          | Green         |  |            |
| SYR1HPL03 | 10 ml      | Concentric          | Orange        |  |            |
| SYR1HPL04 | 20 ml      | Concentric          | Pink          |  | 50         |
| SYR1HPL05 | 50 / 60 ml | Concentric          | Black         |  | 20         |









# Syringes

*Sterilised, Single-use, Luer-Slip*



## Product Features

- Transparent medical polypropylene (PP) barrel.
- Luer-slip tip
- Automated assembly process to prevent personnel contamination.
- Synthetic rubber stopper to reduce allergic reactions.
- Medical grade silicone lubrication to improve sealing ability and glide.
- E.O. Sterilised
- Non-pyrogenic

| Code No.  | Capacity   | Concentr. /Eccentr. | Colour Coding |   | Units /Box |
|-----------|------------|---------------------|---------------|---|------------|
| SYR1HPS01 | 0.5 / 1ml  | Concentric          | Purple        |  | 100        |
| SYR1HPS03 | 2 / 3 ml   | Concentric          | Blue          |  |            |
| SYR1HPS04 | 5 ml       | Concentric          | Green         |  |            |
| SYR1HPS05 | 10 ml      | Concentric          | Orange        |  |            |
| SYR1HPS06 | 20 ml      | Concentric          | Pink          |  | 50         |
| SYR1HPS07 | 50 / 60 ml | Concentric          | Blue          |  | 50         |



# Syringes



*Sterilised, Single-use*



## Product Features – Catheter Syringe

- Transparent medical polypropylene (PP) barrel.
- Catheter tip
- Automated assembly process to prevent personnel contamination.
- Synthetic rubber stopper to reduce allergic reactions.
- Medical grade silicone lubrication to improve sealing ability and glide.
- E.O. Sterilised
- Non-pyrogenic
- 25 Pieces per box


| Code No.  | Capacity | Box /Carton |
|-----------|----------|-------------|
| SYR1CAT01 | 50/60 ml | 16          |

## Product Features – BCG Syringes

- Highly transparent medical polypropylene (PP) barrel
  - Piston retraction lock to prevent re-use
  - E.O. Sterilised
  - Non-pyrogenic
  - 100 pieces per pack
- Manufactured to ISO 8537 requirements

| Code No.  | Capacity | Gauge    | Packs /Carton |
|-----------|----------|----------|---------------|
| SYR1BCG01 | 0,5 ml   | 27G x ½" | 20            |





# Personal, Protective Equipment

*Apparels and Gloves*



Product Features

- Examination gloves
- Powdered or Powder-free
- Finger textured or T-Grip
- Light gauge, Normal or High-risk
- Ambidextrous
- Various cuff lengths
- Single use only
- Packs of 100/box, 10 boxes/case

| Code No.    | Weight | Powder free | Cuff length | Sizes       | Colour Coding |
|-------------|--------|-------------|-------------|-------------|---------------|
| GLV0ELF01-3 | 5,5G   | ✓           | 240mm       | S; M; L     | White         |
| GLV0ELP01-3 | 5G     | ×           | 240mm       | S; M; L     | White         |
| GLV0ELF71-2 | 14G    | ✓           | 400mm       | M; L        | White         |
| GLV0ELF73-6 | 7,5G   | ✓           | 300mm       | S; M; L; XL | White         |
| GLV0EHL01-4 | 17G    | ✓           | 300mm       | S; M; L; XL | Blue          |

# Nitrile gloves

Examination grade

Examtex®



## Product Features

- Examination gloves
- Powdered or Powder-free
- Finger textured or T-Grip
- Light gauge, Normal or High-risk
- Ambidextrous
- Various cuff lengths
- Single use only
- Packs of 100/box, 10 boxes/case

| Code No.     | Weight | Powder free | Cuff length | Sizes           | Colour Coding |
|--------------|--------|-------------|-------------|-----------------|---------------|
| GLV0ENF01-5  | 3,5G   | ✓           | 240mm       | XS; S; M; L; XL | Blue          |
| GLV0ENF71-3  | 21G    | ✓           | 500mm       | M; L; XL        | White         |
| GLV0ENF74-5  | 8,6G   | ✓           | 240mm       | M; L            | Orange        |
| GLV0ENF76-8  | 5G     | ✓           | 240mm       | S; M; L         | Black         |
| GLV0EHN71-72 | 11,4G  | ✓           | 300mm       | M; L            | Blue          |
| GLV0ENF79-80 | 7,7G   | ✓           | 300mm       | M; L            | Black         |
| GLV0ENF81-83 | 5,4G   | ✓           | 240mm       | S; M; L         | Orange        |
| GLV0ENF06-8  | 3,5G   | ✓           | 240mm       | S; M; L         | Violet        |





# Protective wear

Hygienic & Protective



## Product Features – P.E. Gloves


- Smooth film Polyethylene
- Cuff length
- Blue and Transparent
- Ambidextrous
- 0,016cm x 65cm x 110cm
- 100 pieces per pack

| Code No.  | Material | Colour      | Packs /Carton |
|-----------|----------|-------------|---------------|
| GLV0PES02 | LDPE     | Blue        | 100           |
| GLV0PES03 |          | Transparent |               |

## Product Features – P.E. Shoulder Length Gloves

- Light Blue Polyethylene
  - Shoulder length
  - Ambidextrous
- Sterile examination variations available
  - 830mm x 225mm

| Code No.  | Sterile | Units /Pack | Packs /Carton |
|-----------|---------|-------------|---------------|
| GLV0PEL02 | ×       | 50          | 21            |
| GLV1PEL01 | ✓       | 2           | 200           |



# Protective wear

Hygienic & Protective



## Product Features – H.D.P.E. Gloves

- Embossed film High-Density Polyethylene
- Cuff length
- Blue and Transparent
- Ambidextrous
- 100 pieces per pack

| Code No.  | Material | Colour      | Packs /Carton |
|-----------|----------|-------------|---------------|
| GLV0PEH01 | HDPE     | Transparent | 100           |
| GLV0PEH02 |          | Blue        |               |





# Protective wear

Hygienic & Protective



## Product Features – Aprons

- Single use LDPE/HDPE Aprons and LDPE Hang-up Aprons.
- Knee-length with waist ties for secure fit
- 100 pieces per pack for discrete aprons
- 20 pieces per board, 5 boards per pack for Hang-up Aprons
- LDPE aprons in White, Blue, Red, Yellow, Green and Transparent

| Code No.    | Material | Size (cm) | Packs /Carton |
|-------------|----------|-----------|---------------|
| APN0HPE01   | HDPE     | 65 x 110  | 20            |
| APN0LPE01-6 | LDPE     | 65 x 110  | 10            |
| APN0LPH01   | LDPE     | 65 x 110  | 2             |

## Product Features – Smocks

- LDPE Smocks
- Knee-length
- 50 pieces per pack
- White

| Code No.  | Material | Size     | Units /Carton |
|-----------|----------|----------|---------------|
| APN0SMK01 | LDPE     | 80 x 118 | 10            |
| APN0SMK02 | LDPE     | 70 x 118 | 10            |



# Protective wear

Hygienic & Protective



## Product Features – Sleeve Protectors

- Single use LDPE sleeve protectors in premium and regular grades
- 100 pieces per pack, 20 packs per case
- Sleeves in white and blue
- Handmade and Machine-made

| Code No.  | Colour | Manuf. process | Grade   |
|-----------|--------|----------------|---------|
| SLP0PES01 | White  | Machine made   | Regular |
| SLP0PES02 | White  | Machine made   | Premium |
| SLP0PES03 | Blue   | Handmade       | Regular |
| SLP0PES04 | White  | Handmade       | Regular |

## Product Features – Wrist Watch Cover

- Non-woven polypropylene
  - Elasticated
  - White
- 100 pieces per pack

| Code No.  | Colour | Packs /Carton |
|-----------|--------|---------------|
| COV0WTC01 | White  | 10            |





# Protective wear

Hygienic & Protective



## Product Features – Overshoes

- Non-woven polypropylene
- 35 gsm / 25 gsm
- Locally manufactured
- Elasticated for secure fit
- 100 pieces per pack

| Code No.  | Colour | Fabric Weight | Packs /Carton |
|-----------|--------|---------------|---------------|
| COV0SHE02 | Blue   | 25 gsm        | 30            |
| COV0SHE08 | White  | 25 gsm        | 30            |
| COV0SHE10 | Black  | 25 gsm        | 30            |

## Product Features – Boot Covers

- LDPE Single-use boot covers
  - Blue
- Elasticated for Secure fit
- 100 pieces per pack

| Code No.  | Colour | Packs /Carton |
|-----------|--------|---------------|
| COV0SHE06 | Blue   | 10            |





# Protective wear

Hygienic & Protective



## Product Features – Masks

- 3 Ply, Non-woven polypropylene
- Curved side welds to improve seal
- Earloop
- SANS 1886-1:2018 Certified, Medical grade mask
- Type II Bacterial Filtration Efficiency > 98%
- Locally manufactured general use masks. Packed in either pouches of 5 or boxes of 50

| Code No.  | Fastening Method | Colour | Units /Pack | Box /Cart on |
|-----------|------------------|--------|-------------|--------------|
| MSK0NWN02 | Earloop          | Blue   | 50          | 40           |
| MSK0NWN09 | Earloop          | Black  | 50          | 40           |
| MSK0NWN03 | Tie-tape         | Blue   | 50          | 40           |
| MSK0NWN05 | Earloop          | Blue   | 5           | 100          |

## Product Features – Beard Cover

- Non-woven polypropylene
- Ear-loop
- 100 pieces per case


| Code No.  | Size       | Units /Carton |
|-----------|------------|---------------|
| COV0BER01 | 34 x 26 cm | 10            |





# Protective wear

Hygienic & Protective



## Product Features – Scrub Suit

- Non-woven
- Gown & Pant combination
- SMS Grade
- Round collar, elasticated trouser waist
- Individually packed

| Code No.  | Size     | Units /Carton |
|-----------|----------|---------------|
| GWN0SBS01 | Small    | 50            |
| GWN0SBS06 |          | 15            |
| GWN0SBS02 | Medium   | 50            |
| GWN0SBS03 | Large    |               |
| GWN0SBS04 | X-Large  |               |
| GWN0SBS05 | XX-Large |               |

## Product Features – SMS Gowns

- Non-woven polypropylene
  - Hand wipes included
  - SMS grade
- Reinforced and Non-reinforced surgical gowns
  - Waist ties for secure fit
  - Individually packed

| Code No.  | Enforced | Size     | Units /Carton |
|-----------|----------|----------|---------------|
| GWN0SUR01 | ×        | Large    | 20            |
| GWN1RSR01 | ✓        | Medium   |               |
| GWN1RSR02 | ✓        | Large    |               |
| GWN1RSR03 | ✓        | X-Large  |               |
| GWN1RSR04 | ✓        | XX-Large |               |



# Protective wear

Hygienic & Protective



## Product Features – Examination Gown

- Non-woven, spun-bound polypropylene
- Short Sleeve
- Waist tie for secure fit
- Individually packed

| Code No.  | Size    | Packs /Carton |
|-----------|---------|---------------|
| GWN0EXM03 | Large   | 100           |
| GWN0EXM04 | X-Large |               |

## Product Features – Pediatric Examination Gown

- Non-woven, spun-bound polypropylene
- Short Sleeve
- Waist tie for secure fit
- Individually packed

| Code No.  | Material  | Packs /Carton |
|-----------|-----------|---------------|
| GWN0EXM01 | Non-woven | 100           |







# Bonnet cap

Hygienic & Protective



## Product Features

- Non-woven
- Packed in 100 pieces
- Available in Green or White

| Code No.    | Material  | Packs /Carton |
|-------------|-----------|---------------|
| CAP0BON03-4 | Non-woven | 20            |

## Product Features – Mop/Bouffant Caps

- Non-woven
- Elasticated
- Locally manufactured
- Offering sizes in 20 inch, 18 inch and 16 inch
- Packed in 100 pieces
- Available in Black, Blue, White, Yellow, Orange , Red, Pink or Green.

| Code No.      | Material  | Units /Carton |
|---------------|-----------|---------------|
| CAP0MOP01 - 9 | Non-woven | 20            |




A close-up photograph of a hand applying a white, textured bandage to a person's foot. The background is a blurred blue fabric. A blue rectangular overlay with a thin gold border is positioned over the lower part of the image, containing the text 'Wound Care'.

# Wound Care

*Bandages & Dressings*

# Gauze swabs

Centromed



## Product Features

- 100% Bleached cotton, high absorption material.
- Minimal linting design
- E.O. Sterilised and Non-Sterile
- X-ray and plain
- Superior air permeability and softer compared to traditional gauze
- Automated manufacturing process
- Packed in 100 pieces, bulk pack of 200 pieces per pack for 100mm x 100mm 4-ply available

| Code No.     | Size (mm) | Ply | Sterilised | X-Ray | Units /Pack | Units /Carton |
|--------------|-----------|-----|------------|-------|-------------|---------------|
| GZE0W0801    | 50 x 50   | 8   | ×          | ×     | 100         | 40            |
| GZE0W0802    | 75 x 75   | 8   | ×          | ×     |             |               |
| GZE0W0804/05 | 100 x 100 | 8   | ×          | ×     |             |               |
| GZE0W1201    | 100 x 100 | 12  | ×          | ×     |             |               |
| GZE0X0803    | 100 x 100 | 8   | ×          | ✓     | 100         | 40            |
| GZE0X0801    | 50 x 50   | 8   | ×          | ✓     |             |               |
| GZE0X0802    | 75 x 75   | 8   | ×          | ✓     |             |               |
| GZE0X1202    | 100 x 100 | 12  | ×          | ✓     |             |               |
| GZE0X1201    | 75 x 75   | 12  | ×          | ✓     | 100         | 40            |
| GZE0X1601    | 100 x 100 | 16  | ×          | ✓     |             |               |
| GZE0X1602    | 100 x 350 | 16  | ×          | ✓     |             |               |
| GZE0X3201    | 100 x 100 | 32  | ×          | ✓     |             |               |
| GZE1W0803    | 100 x 100 | 8   | ✓          | ×     | 5           | 520           |

# Gauze rolls

Centromed



## Product Features – Gauze Rolls

- Non-sterile
- 100% Bleached cotton, highly absorbent
- Minimal linting & loose end design
- Individually packed

| Code No.  | Size        | Rolls /Carton |
|-----------|-------------|---------------|
| GZE0WRL01 | 90cm x 5m   | 360           |
| GZE0WRL02 | 90cm x 100m | 20            |

## Product Features – Non-woven Gauze Swabs

- Non-woven rayon and polyester synthetic fibre
  - Minimal linting design
  - High absorption
- Superior air permeability and softer compared to traditional gauze
  - Automated manufacturing process
- Packed in 100 pieces, bulk pack of 200 pieces per pack for 100mm x 100mm 4-ply available.

| Code No.    | Size (mm) | Ply | Packs /Carton |
|-------------|-----------|-----|---------------|
| GZE0N0201   | 100 x 100 | 2   | 40            |
| GZE0N0401   | 50 x 50   | 4   | 80            |
| GZE0N0402   | 75 x 75   |     | 40            |
| GZE0N0403/4 | 100 x 100 |     |               |



# Abdominal swabs

Centromed



## Product Features

- 100% Bleached cotton
- Minimal linting & loose end design
- E.O. Sterilised and Non-Sterile
- Single X-ray detectable line
- Single strap design for grip

| Code No.  | Size (mm) | Ply | Sterilised | Units /Pack | Units /Carton |
|-----------|-----------|-----|------------|-------------|---------------|
| GZE0AB401 | 200 x 170 | 4   | ×          | 25          | 32            |
| GZE0AB402 | 350 x 350 | 4   | ×          |             |               |
| GZE0AB403 | 450 x 370 | 4   | ×          |             |               |
| GZE0AB601 | 350 x 350 | 6   | ×          |             |               |
| GZE0AB602 | 450 x 370 | 6   | ×          |             |               |
| GZE1AB601 | 450 x 370 | 6   | ✓          | 5           | 24            |





# Medical Sundries

*Sundry consumables*



# Medical sundries

Sundry consumables



## Product Features – Urine Bags

- A3 (T-type valve) and A4 (Sample valve)
- Type 2 (Double ended with Bottom drain), non-body-worn urine drainage bag.
- Medical grade non-toxic PVC (Polyvinyl Chloride)
- Graduated with transparent sides.
- Non-return valve at inlet tube.
- Tube internal diameters not less than 5mm.

| Code No.  | Class | Capacity(l) | Packs /Carton |
|-----------|-------|-------------|---------------|
| BAG1URN01 | A3    | 2           | 200           |
| BAG1URN02 | A4    | 2           | 200           |

## Product Features – Body Bags

- Meets SA governmental regulations
  - 4 Reinforced grab-handles
    - Name tag pouch
  - Lockable zipper pull
  - U-Shaped opening
- Liquid absorption pad included
  - Individually packed

| Code No.  | Size     |           | Packs /Carton |
|-----------|----------|-----------|---------------|
| BBG0STD08 | XX-Large | 2,4 x 1 m | 12 or 25      |





# Cleaning wipes

Non-woven



## Product Features – Cleaning Wipe Rolls

- All purpose cleaning wipes
- Lint-free non-woven
- Absorbent, washable.
- Perforated, 24 x 50 cm per square.
- Individually packed

| Code No.  | Size (m)   | Rolls /Carton |
|-----------|------------|---------------|
| WPE0CLR02 | 0,24 x 400 | 5             |

## Product Features – Antibacterial Wet-wipes

- Spunlace non-woven fabric
  - Alcohol free
  - Fragrance-free
  - Paraben-free
- Re-sealable opening
  - SGS Certified
- 40 Wipes per pack

| Code No.  | Size (cm) | Packs /Carton |
|-----------|-----------|---------------|
| WPE1WAF01 | 20 x 15   | 24            |





# Medical sundries



## Product Features – Disposable bags

- LDPE smooth
- Transparent red
- 100 pieces per pack

| Code No.  | Size (cm) | Packs /Carton |
|-----------|-----------|---------------|
| BAG0LPE01 | 28 x 40   | 30            |

## Product Features – Infrared Thermometer

- Large LCD screen with backlighting to improve legibility
- Colour-coded backlight to ease identification of abnormal body temperatures
  - Takes AAA batteries
  - CE, FCC and RoHS compliant
  - Individually packed

| Code No.  | Brand  | Packs /Carton |
|-----------|--------|---------------|
| THM0GEN07 | Contec | 50            |



# Medical sundries

Hygienic & Protective



## Product Features – Biohazard Bag

- LDPE smooth
- Transparent red
- 50 Microns
- 50 pieces per pack

| Code No.   | Size (cm) | Packs /Carton |
|------------|-----------|---------------|
| BAG0LPE002 | 50 x 60   | 10            |





## **Central Medical Head Office**

Durban, South Africa

28 Duiker Road, Canelands, 4339

Tel: +27 32 533-9724

P.O Box: 164, Canelands, 4344

**[www.centromed.co.za](http://www.centromed.co.za)**

