



2024

Product Catalogue



Company Profile

About Us

Central Medical (Pty) Ltd was established in 1996 as a local manufacturer of disposable non-woven products, arguably one of the first to introduce automated manufacturing processes for medical supplies into South Africa.

The Company has since diversified its operations to also perform direct importation of additional hygienic disposable products. Over the years the Company has developed a number of internal brands that has accumulated significant reputation for being reliable, consistent and affordable. These brands include LYSA, Centromed, Hygienic & Protective (H&P) and Examtex.

What we can help you with

In addition to providing goods currently listed in our product catalogue, we are able to provide you with medical supplies in your specified branding as well as any medical supplies that would further suite your needs.

What we stand for

Comprehensive understanding of our offerings


Central Medical has manufactured latex gloves, non-woven apparels and bandages locally. In addition, we have long term relationships with our overseas suppliers since the company's inception. The experience held by our team has allowed us to ensure that the products manufactured would meet manufacturing and regulatory standards.

Consistent offerings

We uphold the same manufacturing standards which we implement locally to our overseas long-term suppliers in order to ensure that the quality of our products remain consistent.

Affordable and Competitive

As a manufacturer and a direct importer, Central Medical avoids the additional costs charged by agents. This cost benefit is reflected in our prices, meaning that it is passed onto our partners.


CEO



Contact

Central Medical Branches

Branches



Durban – Head Office

28 Duiker Road, Canelands, 4339

sales@centromed.co.za

T: +27 32 533 9724

C: +27 82 800 6666



Johannesburg – Branch

88 Cranberry Crescent, Unit 6 Acacia Park

sales-jhb@centromed.co.za

C: +27 73 493 2625

C: +27 64 065 3470



2

Products

2024 Offerings



Needles & Syringes

Sterile, Single-use only

Hypodermic Needle

Sterilised, Single-use



Product Features

- Thin-walled needles
- Made of austenite stainless steel
- Smooth surface with light lubricative coating
- Minimal pain upon puncture
- Transparent Luer-lock plastic hub
- Color-coded, medical grade plastic hubs for easy identification
- E.O. Sterilised
- Non-pyrogenic
- Box of 100 pieces

Code No.	Gauge x Length	Ø x Length (mm)	Colour Coding	Units /Carton
NED1HPD01	15G x 1,5"	1,8 x 40	Blue-grey	100
NED1HPD02	18G x 1,5"	1,2 x 40	Pink	
NED1HPD03	19G x 1,5"	1,1 x 40	Cream	
NED1HPD04	20G x 1,5"	0,9 x 40	Yellow	
NED1HPD05	21G x 1,5"	0,8 x 40	Deep-green	
NED1HPD06	22G x 1,25"	0,7 x 32	Black	
NED1HPD07	23G x 1"	0,6 x 25	Deep-blue	
NED1HPD09	25G x 0,625"	0,5 x 16	Orange	
NED1HPD10	26G x 0,625"	0,45 x 16	Brown	
NED1HPD11	29G x 0,33"	0,33 x 8	Red	

Needles



Sterilised, Single-use





Product Features - Blood Lancets

- Twist-off protective cap
- Pressure activated
- Gamma Sterilised
- Single-use only
- 100 pieces per pack

Code No.	Gauge	Colour Coding	Units /Carton
LAN1ECO04	21G * 1,8mm	Green 	20
LAN1ECO02	25G * 1,6mm	Yellow 	
LAN1ECO03	28G * 1,4mm	Orange 	

Product Features – ReadyLance Blood Lancets

- Twist-off protective cap
- Pressure activated
- Gamma Sterilised
- Single-use only
- 100 pieces per pack

Code No.	Gauge	Colour Coding	Units /Carton
LAN1STD01	25G * 1,6mm	Yellow 	20
LAN1STD02	28G * 1,4mm	Orange 	



Needles



Sterilised, Single-use

Product Features – TB Syringes

- Thin-walled long-point needle
- Highly transparent medical polypropylene (PP) barrel
- E.O. Sterilised
- Non-pyrogenic
- Manufactured to ISO 8537 requirements
- 100 pieces per pack

Code No.	Capacity	Gauge	Packs /Carton
SYR1TBA01	1 ml	25G x 5/8"	20



Insulin syringes

Sterilised, Single-use



Product Features

- Double capped to ensure sterility with ultra-thin needles
- No dead-space within syringe cavity to reduce insulin wastage
- Ultra thin, sharp and lubricated needle.
- E.O. sterilised
- Non-pyrogenic
- U-100/U-50
- Manufactured to ISO 8537 requirements
- 100 pieces per pack

Code No.	Capacity	Gauge	Packs /Carton
SYR1INS01	0,3 ml	29G x ¼"	20
SYR1INS02	0,5 ml	29G x ¼"	
SYR1INS03	1,0 ml	29G x ¼"	
SYR1INS04	0,5 ml	27G x ½"	
SYR1INS05	1,0 ml	27G x ½"	
SYR1INS06	1,0 ml	29G x ½"	
SYR1INS07	1,0 ml	26G x ½"	

Syringes

Sterilised, Single-use, Luer-Lock



Product Features

- Transparent medical polypropylene (PP) barrel.
- Luer-lock tip
- Automated assembly process to prevent personnel contamination.
- Synthetic rubber stopper to reduce allergic reactions.
- Medical grade silicone lubrication to improve sealing ability and glide.
- E.O. Sterilised
- Non-pyrogenic

Code No.	Capacity	Concentr. /Eccentr.	Colour Coding	Units /Box
SYR1HPL01	2 / 3 ml	Concentric	Blue ●	100
SYR1HPL02	5 ml	Concentric	Green ●	
SYR1HPL03	10 ml	Concentric	Orange ●	
SYR1HPL04	20 ml	Concentric	Pink ●	50
SYR1HPL05	50 / 60 ml	Concentric	Black ●	20







Syringes

Sterilised, Single-use, Luer-Slip



Product Features

- Transparent medical polypropylene (PP) barrel.
- Luer-slip tip
- Automated assembly process to prevent personnel contamination.
- Synthetic rubber stopper to reduce allergic reactions.
- Medical grade silicone lubrication to improve sealing ability and glide.
- E.O. Sterilised
- Non-pyrogenic

Code No.	Capacity	Concentr. /Eccentr.	Colour Coding	Units /Box
SYR1HPS01	0.5 / 1ml	Concentric	Purple 	100
SYR1HPS03	2 / 3 ml	Concentric	Blue 	
SYR1HPS04	5 ml	Concentric	Green 	
SYR1HPS05	10 ml	Concentric	Orange 	
SYR1HPS06	20 ml	Concentric	Pink 	50
SYR1HPS07	50 / 60 ml	Concentric	Blue 	50

Syringes



Sterilised, Single-use



Product Features – Catheter Syringe

- Transparent medical polypropylene (PP) barrel.
- Catheter tip
- Automated assembly process to prevent personnel contamination.
- Synthetic rubber stopper to reduce allergic reactions.
- Medical grade silicone lubrication to improve sealing ability and glide.
- E.O. Sterilised
- Non-pyrogenic
- 25 Pieces per box


Code No.	Capacity	Box /Carton
SYR1CAT01	50/60 ml	16

Product Features – BCG Syringes

- Highly transparent medical polypropylene (PP) barrel
 - Piston retraction lock to prevent re-use
 - E.O. Sterilised
 - Non-pyrogenic
- 100 pieces per pack
- Manufactured to ISO 8537 requirements

Code No.	Capacity	Gauge	Packs /Carton
SYR1BCG01	0,5 ml	27G x ½"	20





Personal, Protective Equipment

Apparels and Gloves

Latex gloves

Examination grade

Examtex®



Product Features

- Examination gloves
- Powdered or Powder-free
- Finger textured or T-Grip
- Light gauge, Normal or High-risk
- Ambidextrous
- Various cuff lengths
- Single use only
- Packs of 100/box, 10 boxes/case

Code No.	Weight	Powder free	Cuff length	Sizes	Colour Coding
GLV0ELF01-3	5,5G	✓	240mm	S; M; L	White ●
GLV0ELP01-3	5G	×	240mm	S; M; L	White ●
GLV0ELF71-2	14G	✓	400mm	M; L	White ●
GLV0ELF73-6	7,5G	✓	300mm	S; M; L; XL	White ●
GLV0EHL01-4	17G	✓	300mm	S; M; L; XL	Blue ●

Nitrile gloves

Examination grade

Examtex[®]



Product Features

- Examination gloves
- Powdered or Powder-free
- Finger textured or T-Grip
- Light gauge, Normal or High-risk
- Ambidextrous
- Various cuff lengths
- Single use only
- Packs of 100/box, 10 boxes/case

Code No.	Weight	Powder free	Cuff length	Sizes	Colour Coding
GLV0ENF01-5	3,5G	✓	240mm	XS; S; M; L; XL	Blue 
GLV0ENF71-3	21G	✓	500mm	M; L; XL	White 
GLV0ENF74-5	8,6G	✓	240mm	M; L	Orange 
GLV0ENF76-8	5G	✓	240mm	S; M; L	Black 
GLV0EHN71-72	11,4G	✓	300mm	M; L	Blue 
GLV0ENF79-80	7,7G	✓	300mm	M; L	Black 
GLV0ENF81-83	5,4G	✓	240mm	S; M; L	Orange 
GLV0ENF06-8	3,5G	✓	240mm	S; M; L	Violet 

Protective wear

Hygienic & Protective



Product Features – P.E. Gloves

- Smooth film Polyethylene
- Cuff length
- Blue and Transparent
- Ambidextrous
- 0,016cm x 65cm x 110cm
- 100 pieces per pack

Code No.	Material	Colour	Packs /Carton
GLV0PES02	LDPE	Blue	100
GLV0PES03		Transparent	

Product Features – P.E. Shoulder Length Gloves

- Light Blue Polyethylene
 - Shoulder length
 - Ambidextrous
- Sterile examination variations available
 - 830mm x 225mm

Code No.	Sterile	Units /Pack	Packs /Carton
GLV0PEL02	×	50	21
GLV1PEL01	✓	2	200



Protective wear

Hygienic & Protective



Product Features – H.D.P.E. Gloves

- Embossed film High-Density Polyethylene
- Cuff length
- Blue and Transparent
- Ambidextrous
- 100 pieces per pack

Code No.	Material	Colour	Packs /Carton
GLV0PEH01	HDPE	Transparent	100
GLV0PEH02		Blue	



Protective wear

Hygienic & Protective



Product Features – Aprons

- Single use LDPE/HDPE Aprons and LDPE Hang-up Aprons.
- Knee-length with waist ties for secure fit
- 100 pieces per pack for discrete aprons
- 20 pieces per board, 5 boards per pack for Hang-up Aprons
- LDPE aprons in White, Blue, Red, Yellow, Green and Transparent

Code No.	Material	Size (cm)	Packs /Carton
APN0HPE01	HDPE	65 x 110	20
APN0LPE01-6	LDPE	65 x 110	10
APN0LPH01	LDPE	65 x 110	2

Product Features – Smocks

- LDPE Smocks
- Knee-length
- 50 pieces per pack
 - White

Code No.	Material	Size	Units /Carton
APN0SMK01	LDPE	80 x 118	10
APN0SMK02	LDPE	70 x 118	10



Protective wear

Hygienic & Protective



Product Features – Sleeve Protectors

- Single use LDPE sleeve protectors in premium and regular grades
- 100 pieces per pack, 20 packs per case
- Sleeves in white and blue
- Handmade and Machine-made

Code No.	Colour	Manuf. process	Grade
SLP0PES01	White	Machine made	Regular
SLP0PES02	White	Machine made	Premium
SLP0PES03	Blue	Handmade	Regular
SLP0PES04	White	Handmade	Regular

Product Features – Wrist Watch Cover

- Non-woven polypropylene
 - Elasticated
 - White
- 100 pieces per pack

Code No.	Colour	Packs /Carton
COV0WTC01	White	10



Protective wear

Hygienic & Protective



Product Features – Overshoes

- Non-woven polypropylene
- 35 gsm/25 gsm
- Locally manufactured
- Elasticated for secure fit
- 100 pieces per pack

Code No.	Colour	Fabric Weight	Packs /Carton
COV0SHE02	Blue	25 gsm	30
COV0SHE08	White	25 gsm	30
COV0SHE10	Black	25 gsm	30

Product Features – Boot Covers

- LDPE Single-use boot covers
 - Blue
- Elasticated for Secure fit
 - 100 pieces per pack

Code No.	Colour	Packs /Carton
COV0SHE06	Blue	10



Protective wear

Hygienic & Protective



Product Features – Masks

- 3 Ply, Non-woven polypropylene
- Curved side welds to improve seal
- Earloop
- SANS 1886-1:2018 Certified, Medical grade mask
- Type II Bacterial Filtration Efficiency > 98%
- Locally manufactured general use masks. Packed in either pouches of 5 or boxes of 50

Code No.	Fastening Method	Colour	Units /Pack	Box /Cart on
MSK0NWN02	Earloop	Blue	50	40
MSK0NWN09	Earloop	Black	50	40
MSK0NWN03	Tie-tape	Blue	50	40
MSK0NWN05	Earloop	Blue	5	100

Product Features – Beard Cover

- Non-woven polypropylene
- Ear-loop
- 100 pieces per case

Code No.	Size	Units /Carton
COV0BER01	34 x 26 cm	10



Protective wear

Hygienic & Protective



Product Features – Scrub Suit

- Non-woven
- Gown & Pant combination
- SMS Grade
- Round collar, elasticated trouser waist
- Individually packed

Code No.	Size	Units /Carton
GWN0SBS01	Small	50
GWN0SBS06		15
GWN0SBS02	Medium	50
GWN0SBS03	Large	
GWN0SBS04	X-Large	
GWN0SBS05	XX-Large	

Product Features – SMS Gowns

- Non-woven polypropylene
 - Hand wipes included
 - SMS grade
- Reinforced and Non-reinforced surgical gowns
 - Waist ties for secure fit
 - Individually packed

Code No.	Enforced	Size	Units /Carton
GWN0SUR01	×	Large	20
GWN1RSR01	✓	Medium	
GWN1RSR02	✓	Large	
GWN1RSR03	✓	X-Large	
GWN1RSR04	✓	XX-Large	



Protective wear

Hygienic & Protective



Product Features – Examination Gown

- Non-woven, spun-bound polypropylene
- Short Sleeve
- Waist tie for secure fit
- Individually packed

Code No.	Size	Packs /Carton
GWN0EXM03	Large	100
GWN0EXM04	X-Large	

Product Features – Pediatric Examination Gown

- Non-woven, spun-bound polypropylene
 - Short Sleeve
- Waist tie for secure fit
- Individually packed

Code No.	Material	Packs /Carton
GWN0EXM01	Non-woven	100



Bonnet cap

Hygienic & Protective



Product Features

- Non-woven
- Packed in 100 pieces
- Available in Green or White

Code No.	Material	Packs /Carton
CAP0BON03-4	Non-woven	20

Product Features – Mop/Bouffant Caps

- Non-woven
- Elasticated
- Locally manufactured
- Offering sizes in 20 inch, 18 inch and 16 inch
- Packed in 100 pieces
- Available in Black, Blue, White, Yellow, Orange, Red, Pink or Green.

Code No.	Material	Units /Carton
CAP0MOP01 - 9	Non-woven	20



A close-up photograph of a person's hand applying a white, textured bandage to another person's arm. The background is a blurred blue fabric. A blue rectangular overlay is positioned in the lower-left quadrant, containing the text 'Wound Care'.

Wound Care

Bandages & Dressings

Gauze swabs

Centromed



Product Features

- 100% Bleached cotton, high absorption material.
- Minimal linting design
- E.O. Sterilised and Non-Sterile
- X-ray and plain
- Superior air permeability and softer compared to traditional gauze
- Automated manufacturing process
- Packed in 100 pieces, bulk pack of 200 pieces per pack for 100mm x 100mm 4-ply available

Code No.	Size (mm)	Ply	Sterilised	X-Ray	Units /Pack	Units /Carton
GZE0W0801	50 x 50	8	×	×	25	32
GZE0W0802	75 x 75	8	×	×		
GZE0W0803	100 x 100	8	×	×		
GZE0W1201	100 x 100	12	×	×		
GZE0X0803	100 x 100	8	×	✓		
GZE0X0801	50 x 50	8	×	✓		
GZE0X0802	75 x 75	8	×	✓		
GZE0X1202	100 x 100	12	×	✓		
GZE0X1201	75 x 75	12	×	✓		
GZE0X1601	100 x 100	16	×	✓		
GZE0X1602	100 x 350	16	×	✓		
GZE0X3201	100 x 100	32	×	✓		
GZE1W0803	100 x 100	8	✓	×	5	24

Gauze rolls

Centromed



Product Features – Gauze Rolls

- Non-sterile
- 100% Bleached cotton, highly absorbent
- Minimal linting & loose end design
- Individually packed

Code No.	Size	Rolls /Carton
GZE0WRL01	90cm x 5m	360
GZE0WRL02	90cm x 100m	20

Product Features – Non-woven Gauze Swabs

- Non-woven rayon and polyester synthetic fibre
 - Minimal linting design
 - High absorption
- Superior air permeability and softer compared to traditional gauze
 - Automated manufacturing process
- Packed in 100 pieces, bulk pack of 200 pieces per pack for 100mm x 100mm 4-ply available.

Code No.	Size (mm)	Ply	Packs /Carton
GZE0N0201	100 x 100	2	40
GZE0N0401	50 x 50	4	80
GZE0N0402	75 x 75		40
GZE0N0403/4	100 x 100		



Abdominal swabs

Centromed



Product Features

- 100% Bleached cotton
- Minimal linting & loose end design
- E.O. Sterilised and Non-Sterile
- Single X-ray detectable line
- Single strap design for grip

Code No.	Size (mm)	Ply	Sterilised	Units /Pack	Units /Carton
GZE0AB401	200 x 170	4	×	25	32
GZE0AB402	350 x 350	4	×		
GZE0AB403	450 x 370	4	×		
GZE0AB601	350 x 350	6	×		
GZE0AB602	450 x 370	6	×		
GZE1AB601	450 x 370	6	✓	5	24



Medical Sundries

Sundry consumables

Medical sundries



Sundry consumables



Product Features – Urine Bags

- A3 (T-type valve) and A4 (Sample valve)
- Type 2 (Double ended with Bottom drain), non-body-worn urine drainage bag.
- Medical grade non-toxic PVC (Polyvinyl Chloride)
- Graduated with transparent sides.
- Non-return valve at inlet tube.
- Tube internal diameters not less than 5mm.

Code No.	Class	Capacity(l)	Packs /Carton
BAG1URN01	A3	2	200
BAG1URN02	A4	2	200

Product Features – Body Bags

- Meets SA governmental regulations
 - 4 Reinforced grab-handles
 - Name tag pouch
 - Lockable zipper pull
 - U-Shaped opening
- Liquid absorption pad included
 - Individually packed

Code No.	Size		Packs /Carton
BBG0STD08	XX-Large	2,4 x 1 m	12 or 25



Cleaning wipes



Non-woven



Product Features – Cleaning Wipe Rolls

- All purpose cleaning wipes
- Lint-free non-woven
- Absorbent, washable.
- Perforated, 24 x 50 cm per square.
- Individually packed

Code No.	Size (m)	Rolls /Carton
WPE0CLR02	0,24 x 400	5

Product Features – Antibacterial Wet-wipes

- Spunlace non-woven fabric
 - Alcohol free
 - Fragrance-free
 - Paraben-free
- Re-sealable opening
 - SGS Certified
- 40 Wipes per pack

Code No.	Size (cm)	Packs /Carton
WPE1WAF01	20 x 15	24



Medical sundries



Product Features – Disposable bags

- LDPE smooth
- Transparent red
- 100 pieces per pack

Code No.	Size (cm)	Packs /Carton
BAG0LPE01	28 x 40	30

Product Features – Infrared Thermometer

- Large LCD screen with backlighting to improve legibility
- Colour-coded backlight to ease identification of abnormal body temperatures
 - Takes AAA batteries
 - CE, FCC and RoHS compliant
 - Individually packed

Code No.	Brand	Packs /Carton
THM0GEN07	Contec	50



Medical sundries

Hygienic & Protective



Product Features – Biohazard Bag

- LDPE smooth
- Transparent red
- 50 Microns
- 50 pieces per pack

Code No.	Size (cm)	Packs /Carton
BAG0LPE002	50 x 60	10





Central Medical Head Office

Durban, South Africa
28 Duiker Road, Canelands, 4339

Tel: +27 32 533-9724
P.O Box: 164, Canelands, 4344
www.centromed.co.za

