



Product specification

Examination Gloves- NitriPlus

Functional Benefits & Product Characteristics

- Protection from unwanted or dangerous substance.
- Beaded cuff provides easy wear & deters roll-back.
- Softness provides superior comfort and natural fit
- Ambidextrous and straight fingers.
- 240 mm cuff
- Latex free and powder free

Product information

Type:	Examination NitriPlus
Sterility:	Non-sterile
Material:	Blended nitrile (not 100% nitrile)
Design & Features:	Ambidextrous; smooth or fully textured surface; beaded cuff.
Colour:	Blue
Powder:	Powder free
Storage:	Product must be stored in a dry area. Product must be kept away from extreme temperatures and sources of heat
Shelf-life:	5 years from the date of manufacturing
Product Code:	GLV-0-EBF-02


Physical Dimensions

Size M weight (g):	5g
Fingertip (mm):	0.1 +/- 0.03
Palm (mm):	0.09 +/- 0.02
Cuff (mm):	0.07 +/- 0.02
Elongation (%) before & after aging:	450/350
Tensile strength:	14



Packaging information

- 100 pcs per box, 10 boxes per case



*mock-ups is for illustration purposes. The final product may vary from picture.